


CAPTOR SERIES

CAPTOR
Thrombectomy Device

SUPSELEK
Microcatheter

extraflex
Distal Access Catheter

FULLBLOCK
Balloon Guiding Catheter



ORDERING INFORMATION
FULLBLOCK Balloon Guiding Catheter

Model	Effective Length	Outer Diameter	Inner Diameter	Compatibility
SM*BGC-S80	80±5cm	8F(2.6±0.5mm)	0.085inch	Compatible with 6f intermediate catheter and 8f/9f introducer sheath
SM*BGC-S95	95±5cm	8F(2.6±0.5mm)		
SM*BGC-S+80	80±5cm	9F(3±0.5mm)	0.088inch	Compatible with 6f intermediate catheter and 9f/10f introducer sheath
SM*BGC-S+95	95±5cm	9F(3±0.5mm)		



This brochure is only delivered to specific person for his or her own reference, is not to be distributed widely as an advertisement. Any information contained in this brochure may be changed without prior notice.

SHANGHAI HEARTCARE MEDICAL TECHNOLOGY CO., LTD.
Building 38, No.356 Zhengbo Road, Lingang New Area,
China (Shanghai) Pilot Free Trade Zone
+86-21-58975055
+86-21-58975056
www.heartcare.com.cn

CAPTOR SERIES



Balloon Guiding Catheter




ADAPTABLE AND
TRUSTWORTHY




Good Compliance

The balloon is made of super-compliant material which allows better apposition, reduces vessel injury and has good occlusion effect even at bifurcation.



Compliance Test 1

Red is inflated balloon and blue is simulated blood




Compliance Test 2

Occlusion effect of balloon guide catheter in irregular models


High Inflation Volume

The balloon remains intact and safe when inflated to 4 times of the standard inflation volume.




Safety Guaranteed

Uniquely laser treated atraumatic tip



Good Compatibility

Over 6F inner diameter



Efficient Occlusion

Super compliance balloon material