

# SC HONKYTONK™

## Technical Information

Balloon Material	PEBAX
Shaft Diameters	2.0F proximal, 2.5F distal
Usable Length	145 cm
Minimum Guiding Catheter	5F (0.056")
Maximum Guide Wire	0.014"

## Order Information

Diameter (mm)	Length (mm)						
	9	12	15	20	25	30	35
1.25	SCB-1.25-09-RX	/	SCB-1.25-15-RX	SCB-1.25-20-RX	/	/	/
1.50	SCB-1.50-09-RX	/	SCB-1.50-15-RX	SCB-1.50-20-RX	/	/	/
2.00	SCB-2.00-09-RX	/	SCB-2.00-15-RX	SCB-2.00-20-RX	SCB-2.00-25-RX	/	/
2.25	SCB-2.25-09-RX	/	SCB-2.25-15-RX	SCB-2.25-20-RX	SCB-2.25-25-RX	SCB-2.25-30-RX	SCB-2.25-35-RX
2.50	SCB-2.50-09-RX	SCB-2.50-12-RX	SCB-2.50-15-RX	SCB-2.50-20-RX	SCB-2.50-25-RX	SCB-2.50-30-RX	SCB-2.50-35-RX
2.75	SCB-2.75-09-RX	SCB-2.75-12-RX	SCB-2.75-15-RX	SCB-2.75-20-RX	SCB-2.75-25-RX	SCB-2.75-30-RX	SCB-2.75-35-RX
3.00	SCB-3.00-09-RX	SCB-3.00-12-RX	SCB-3.00-15-RX	SCB-3.00-20-RX	SCB-3.00-25-RX	SCB-3.00-30-RX	SCB-3.00-35-RX
3.50	SCB-3.50-09-RX	SCB-3.50-12-RX	SCB-3.50-15-RX	SCB-3.50-20-RX	SCB-3.50-25-RX	SCB-3.50-30-RX	SCB-3.50-35-RX
4.00	SCB-4.00-09-RX	SCB-4.00-12-RX	SCB-4.00-15-RX	SCB-4.00-20-RX	SCB-4.00-25-RX	SCB-4.00-30-RX	SCB-4.00-35-RX

## Compliance Chart

Pressure (atm)	Diameter (mm)								
	1.25	1.50	2.00	2.25	2.50	2.75	3.00	3.50	4.00
6	1.25	1.50	2.00	2.25	2.50	2.75	3.00	3.50	4.00
12	1.40	1.64	2.21	2.54	2.74	3.09	3.28	3.85	4.33
14	1.45	1.69	2.28	2.62	2.81	3.17	3.36	3.96	/

NP RBP

© 2022 SINOMED - All rights reserved

SINOMED and NC Rockstar are registered trademarks of Sino Medical Sciences Technology Inc. in USA.

**Sino Medical Sciences Technology Inc.**

2nd Floor, TEDA Biopharm Res, Building B, #5 4th St, TEDA, Tianjin, China

Z: 300457


T: + 86 022 5986 2900

F: + 86 022 5986 2904



# SC HONKYTONK™

## PTCA Balloon Dilatation Catheter



- ✓ Ultimate Crossability
- ✓ Optimized Trackability
- ✓ Full Size Kissing Compatible



## Outstanding Crossing and Tip Entry Profile (Φ 3.0 mm)

0.0284



SINO MED  
Innovation for health

0.0309



Boston Scientific

0.0315




TERUMO

0.0318




Boston Scientific



98%



HYDROPHILIC COATING

Reduces friction by 98% in the coated segment



## Full Size 6F kissing Compatibility


Full size of SC Honkytonk could use more than one balloon in a guiding catheter



## Optimized Trackability

Triple-layer inner tubing and new materials in distal and proximal outer tubes allows for better trackability after multiple uses.

Trackability, for balloon size 2.50 mm x 15 mm



SINO MED  
Innovation for health

Boston Scientific

TERUMO

Boston Scientific