


CAPTOR SERIES

CAPTOR
Thrombectomy Device

SUPSELEK
Microcatheter

extraflex
Distal Access Catheter

FULLBLOCK
Balloon Guiding Catheter



ORDERING INFORMATION
SUPSELEK Micro Catheter

Model	Specification	Effective Length (cm)	Outer Diameter (Fr)	Inner Diameter (inch)
21	SM*MC-21110	110	2.5	0.021
	SM*MC-21130	130	2.5	0.021
	SM*MC-21153	153	2.5	0.021
27	SM*MC-27110	110	2.9	0.027
	SM*MC-27130	130	2.9	0.027
	SM*MC-27145	145	2.9	0.027
	SM*MC-27153	153	2.9	0.027




HeartCare
COMPLETE SOLUTION FOR STROKE
CHERISH HEALTH OF LIFE WITH HEARTCARE

This brochure is only delivered to specific person for his or her own reference, is not to be distributed widely as an advertisement. Any information contained in this brochure may be changed without prior notice.

SHANGHAI HEARTCARE MEDICAL TECHNOLOGY CO., LTD.
Building 38, No.356 Zhengbo Road, Lingang New Area,
China (Shanghai) Pilot Free Trade Zone
+86-21-58975055
+86-21-58975056
www.heartcare.com.cn

CAPTOR SERIES

SUPSELEK
Micro Catheter





HeartCare
COMPLETE SOLUTION FOR STROKE

Flexibility
Performance
Test

The tested maximum yield force for 2mm absolute axial downward shift on a universal tester, with a testing rate of 2mm/min.

© Comparison result: HeartCare SUPSELEK micro catheter had better flexibility than comparison group.



Effective length of 153 cm